

# AlphaTrak<sup>3</sup>

## Demonstrated accuracy and validity for canine, feline and equine results<sup>1,2</sup>

In two prospective clinical studies, the accuracy and clinical relevance of the AlphaTrak<sup>3</sup> glucometer were assessed in comparison with the gold standard reference analyzer: Beckman Coulter<sup>®</sup> Chemistry Analyzer.



### Why use a veterinary glucometer?

The physiological differences in glucose metabolism between humans and veterinary species can lead to inconsistencies in results when using human-specific glucometers.<sup>3</sup>

### Study Design

AlphaTrak 3 comparison with Beckman Coulter Reference Analyzer



### Accuracy Criteria

Studies followed International Organization for Standardization (ISO) guidelines.

AlphaTrak 3 meets the same high requirements for blood-glucose monitoring accuracy criteria used for human glucometers (ISO 15197:2013).

**95%**

of all results must fall within  $\pm 15$  mg/dL (at concentrations  $< 100$  mg/dL) or  $\pm 15\%$  (at concentrations  $\geq 100$  mg/dL) of reference method

**99%**

of blood glucose readings must fall within zones A and B of the consensus error grid, which measures the clinical significance of accuracy results



### Results

Accuracy of AlphaTrak 3 was comparable to Beckman Coulter Reference Analyzer.<sup>1,2</sup>

### AlphaTrak 3 Features & Benefits



Accurate<sup>1,2</sup> and reliable



Small sample size (0.3  $\mu$ L)



Wide testing range (25–750 mg/dL)



Validated for dogs, cats and horses



Bluetooth<sup>®</sup> connected




Mobile App tracks patient glucose levels and activities

# AlphaTrak<sup>®</sup> 3 Clinical Relevance<sup>1,2</sup>

AlphaTrak 3 demonstrated **clinically accurate measurements** across all ranges in the Consensus Error Grid (CEG).



Canine Consensus Error Grid (CEG)




**OBJECTIVE**  
99%

**RESULT**  
**100%**



Feline Consensus Error Grid (CEG)




**OBJECTIVE**  
99%

**RESULT**  
**100%**



Equine Consensus Error Grid (CEG)<sup>†</sup>



**OBJECTIVE**  
99%\*


**RESULT**  
**100%**

■ **ZONE A:** Accurate ■ **ZONE B:** Minor error (no clinical impact) ■ **ZONE C:** Over-correction ■ **ZONE D:** Failure to detect ■ **ZONE E:** Erroneous

AlphaTrak 3 consistently provided **accurate blood glucose readings**, meeting ISO standards for clinical use.



Canine Accuracy (Bias) Plot




**OBJECTIVE**  
95%<sup>†</sup>

**RESULT**  
**95.1%**



Feline Accuracy (Bias) Plot




**OBJECTIVE**  
95%<sup>†</sup>

**RESULT**  
**96.2%**



Equine Accuracy (Bias) Plot



**OBJECTIVE**  
95%<sup>†</sup>

**RESULT**  
**98.9%**

## Summary<sup>1,2</sup>

- ✓ AlphaTrak 3 demonstrated excellent performance when compared to reference results in dogs, cats and horses.
- ✓ AlphaTrak 3 met the highest human ISO accuracy criteria and study protocol criteria for dogs, cats and horses and was comparable to the Beckman Coulter analyzer.
- ✓ 95.1% of canine results, 96.2% of feline results and 98.9% of equine results fell within the ISO accuracy (bias) threshold.
- ✓ 100% of canine, feline and equine results were within zones A and B of the consensus error grid.



**Request a demo of the AlphaTrak 3 today**

\* Equine samples were laboratory manipulated to replicate entire range

<sup>†</sup> Of samples within accuracy threshold

**References:** 1. Data on file. Study No. D876R-US-21-040, 2022, Zoetis Inc.

2. Data on file. Study No. D850R-US-22-005, 2022, Zoetis Inc. 3. Dobromylskyj MJ, Sparkes AH. Assessing portable blood glucose meters for clinical use in cats in the United Kingdom. Vet Rec. 2010 Sep 18;167(12):438-42. PMID: 20852246.

**LOOK DEEPER**

Learn more at [zoetisdiagnostics.com/point-of-care/handhelds/alphatrak](https://zoetisdiagnostics.com/point-of-care/handhelds/alphatrak)

All trademarks are the property of Zoetis Services LLC or a related company or a licensor unless otherwise noted. Bluetooth is a registered trademark of the Bluetooth Special Interest Group. Beckman Coulter is a registered trademark of Beckman Coulter, Inc. © 2024 Zoetis Services LLC. All rights reserved. ATK-00421

**zoetis**