



## Accurate blood glucose monitoring for optimal care<sup>1</sup>

Feel confident partnering with your clients to give their diabetic pets the best possible care. AlphaTrak 3 makes testing and monitoring blood glucose levels easier for everyone.

## AlphaTRAK 2 VS AlphaTrak 3

### Glucose Monitoring System Comparison



	AlphaTrak 2	AlphaTrak 3
Minimum blood volume	0.3 µL	0.3 µL
Blood glucose range	20-750 mg/dL	20-750 mg/dL
Hematocrit range	15-65%	15-65%
Strip turns meter on	Yes	Yes
Time to results	10 seconds	5 seconds
Results stored in meter	250	500
Bluetooth enabled	No	Yes
Web-based App for clinic use	No	Yes
Owner can log pet activity data	No	Yes
Data sharing	No	Yes
FUSE/PIMS/ZoetisDx Enabled	No	Yes
Meets ISO Guidelines*	No <sup>2</sup>	Yes <sup>1</sup>
Animal use	Canine, Feline	Canine, Feline

\* International Organization for Standardization 15197:2013.

## Improve patient outcomes

See how AlphaTrak 3 can help you partner with clients to increase the quality of care for your patients.



1. Data on file, Study Report No. D876R-US-21-040, Zoetis 2022.

2. Wolfenden, G., James, F.E., Hung, L.H.T., Bruce, M. and Thompson, M. (2022), Comparative accuracy of two veterinary-calibrated point-of-care glucometres for measurement of blood glucose concentration in dogs. J Small Anim Pract, 63: 512-519