

Anesthetic Machine Maintenance

Daily Machine Maintenance Procedure

Week of:	Sunday	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday
Primary O ₂ Source							
Back up O ₂ Source							
Oxygen alarm working (if present)							
Flowmeters working							
Vaporizer attached/full							
Leak test performed							
Monitoring equipment functioning							
Active Scavenge system working							
Charcoal Waste Anesthetic Absorber weight							

Changed	CO ₂ Absorber Log (checked box = 1 hour used)				
Date:					
Initials:					

Anesthetic Machine Service Log

Date Serviced	Company/Initials